

**Innovation!**


**Articulating cup design  
for a superior polish**



**new!**  
**ProFlare®**

Articulating Disposable  
Prophy Angle

The new ProFlare® articulating disposable prophy angle provides a better way to polish. Its unique latex-free articulating cup\* pivots and flares adapting to the contour of each tooth surface with less wrist motion for better ergonomics, while its outer ridges reach interproximally. The efficient gear mechanism provides smooth acceleration and optimal jam-free performance. ProFlare offers you superior quality at an affordable price.



**Articulating cup design!**



Unique pivot  
& flare design\*

- Maximum surface adaptation
- Latex-free

\*Patent pending



Contours to  
each tooth  
*More efficient  
polishing*



Less wrist  
motion  
*Better  
ergonomics*



Affordably  
priced  
*Superior  
quality*

**premier®** Inspired solutions  
for daily dentistry®

Premier® Dental Products Company • 888-670-6100 • premierdentalco.com

©2024 Premier Dental | All Rights Reserved

# ProFlare®

Articulating Disposable Propphy Angle



## Innovative articulating cup




**Unique pivot & flare design\***  
• Maximum surface adaptation  
• Latex-free  
*\*Patent pending*



**Low profile**  
• Flares & flexes easily  
• Cleans beneath gingiva



**Contours to each tooth**  
More efficient polishing



**Outer ridges**  
• Interproximal stain removal  
• Microgrooves & ribs for paste retention

## Efficient gear mechanism



**Smooth acceleration**  
Eliminates chatter  
  
**Jam-free**  
Optimal performance



Competitors

ProFlare®

**Compact head & neck**  
Improved vision  
& access




**ProFlare®**

5500150 Bag of 125


Up to  
**← 10% →**  
**greater flaring**

**ProFlare® exhibits up to 10% greater flaring capacity**


Resulting in better stain removal & polish\*




ProFlare®



Brand A



Brand B



Brand C

Actual images taken at 150g of pressure.

\*R. P. Christensen and V. W. Bangarter, Determination of RPM, Time, and Load Used in Oral Prophylaxis Polishing in Vivo; J Dent Res 63(12):1376-1382, December, 1984 / 0224042 Rev0

**premier®** Inspired solutions for daily dentistry®



Premier® Dental Products Company • 1710 Romano Drive, Plymouth Meeting, PA 19462, USA

Phone 888-670-6100 • 610-239-6000 • Fax 610-239-6171 • premierdentalco.com • E-mail: info@premierdentalco.com

Lit# 8625561 0524051 Rev2 PDF