

Introducing!

TSZtech™ ZX Diamond Burs with LTech™ Technology

30% FASTER CUTTING¹

NEW!



TSZTECH™

For Removing & Adjusting Zirconia Restorations

TSZtech™ Diamond Burs are specifically crafted for adjusting and removing zirconia restorations, setting a new standard in precision and efficiency. Now, with the introduction of LTech™ Technology, TSZtech™ ZX is revolutionizing zirconia removal.



30% faster cutting¹,
longer lasting
than other
zirconia burs



Consistent diamond
layer height for
maximum cutting
efficiency



Reduced vibration
for a smoother,
more controlled
operation

Available through your authorized dealer.

Distributed by
premier® Inspired solutions
for daily dentistry®

Premier® Dental Products Company • 888-670-6100 • premierdentalco.com

©2025 Premier Dental | All Rights Reserved


TSZTECH™

Removing & Adjusting Zirconia Restorations


TSZtech™ ZX Diamond Burs with LTech™ Technology

Zirconia removal has long been a challenge, but we're changing that! LTech™ Technology offers a breakthrough in zirconia removal, featuring an advanced contact layer design that prevents diamond dislodgement, ensuring sustained cutting performance.

Paired with tough, engineered crystals that remain sharp and intact even when cutting ultra-hard materials, the TSZtech™ ZX with LTech™ Technology streamlines zirconia removal procedures, saving valuable time and enhancing both patient and dentist experiences.



Proprietary LTech™ Technology
Promotes cleaving of higher points and leveling of the contact edge, without releasing the entire crystals.



Traditional Electroplated Technology
The highest crystals are prone to rapid plucking.

The full line of TSZtech™ Diamond Burs optimize zirconia procedures through advanced grit selection, allowing professionals to achieve desired results with ease.

ZX		2013680		2013681	2013682	2013683	2013684		2013685	2013686
ZM	2015670	2015671	2015672	2015673	2015674	2015675	2015676	2015677	2015678	2015679
ZF	2015660	2015661	2015662	2015663	2015664	2015665	2015666	2015667	2015668	2015669



770.8 Round-End Taper Section & Remove

ZX Featuring LTech™ Proprietary Technology, this bur is designed for fast, efficient zirconia removal, reducing chair time and improving the experience for both patient and dentist.

ZM Ideal for efficient removal and moderate amounts of adjustments, preparing surfaces for the ZF bur.²

ZF Perfect for small, precise final adjustments of zirconia surfaces, before final zirconia polishing.²

TSZtech ZF and ZM diamond burs are available in 5-packs.

TSZtech ZX diamond burs are available in 3-packs.

Distributed by
premier Inspired solutions for daily dentistry®



Premier® Dental Products Company • 1710 Romano Drive, Plymouth Meeting, PA 19462, USA
Phone 888-670-6100 • 610-239-6000 • Fax 610-239-6171 • premierdentalco.com • E-mail: info@premierdentalco.com

Two Striper® and P.B.S.® are registered trademarks and TSZtech™ and LTech™ are trademarks of the manufacturer, Abrasive Technology, LLC. / 1. ZX Grit - Data on file. / 2. Features proprietary P.B.S.® brazing process.

Lit# 8621311 1024058 Rev0 NKP5M