



PremierAir® Barnhart 5/6 Ti • 9061460  
Versatile design for supragingival  
& subgingival cleaning



PremierAir® NV4 Ti • 9061461  
Ideal for interproximal access



PremierAir® Langer 3/4 Ti • 9061462  
Specialized for molar and premolar implant abutments



PremierAir® Neb 128B / Langer 5 Mini Ti • 9061463  
• N128 sickle scaler for narrow spaces  
• Langer 5 universal anterior curette for upper & lower anterior implants

## PremierAir® Titanium Implant Scalers

*For Complete Implant Maintenance*

Premier® Dental Products presents the PremierAir® line of Titanium Implant Scalers, expertly crafted for comprehensive implant maintenance. These medical-grade, high-quality instruments are designed specifically for debridement and biofilm removal with supra and subgingival capabilities.

*Available through your authorized dealer.*



Scan to view full line of PremierAir® instruments!



**Optimal working edges**  
To remove plaque, calculus, and biofilm safely & effectively



**High quality**  
Medical grade titanium alloy working ends; matching implant material for safe, non-scratching maintenance of implant abutments and threads



**Ergonomic excellence**  
Lightweight, tactile resin PremierAir handles. Cool gray color for easy recognition




**Extended longevity**  
Working edges can be resharpened to help maximize lifespan

## Premier® Implant Scalars


Engineered for reliable implant maintenance

- ✓ Fiber reinforced graphite – *Will not scratch implant structures*
- ✓ Precisely engineered tips – *Thin and sharp with well-defined edges, effective in removing plaque and calculus*
- ✓ Reusable – *Autoclavable (273°F / 134°C)*



## Premier® PerioWise® "The Friendly Probe"®

- ✓ Slightly flexible for greater patient comfort
- ✓ Clearly visible when using intra-oral camera
- ✓ Safe around implants
- ✓ UNC 12 provides precision measuring of pocket depths with markings indicating 1mm increments up to 12mm



## Premier® Implant Cement™ Plus

Non-Eugenol, Temporary Resin Cement for Implant Retained Crowns

- ✓ Easy to use – *no etch, primers or adhesive required*
- ✓ Unsurpassed radiopacity – *excellent visibility on radiographs*
- ✓ Secure retention – *for long-lasting hold*
- ✓ Retrievable – *without damage to the crown, abutment or implant*
- ✓ Gel stage in just 2 minutes – *excess easy to remove*



**Standard Pack** - 1-5ml Automix Syringe, 10 Mix Tips  
 3001550 Shade Original (Pink) 3001551 Shade White

**Value Pack** - 3-5ml Automix Syringes, 25 Mix Tips  
 3001555 Shade Original (Pink) 3001556 Shade White

**Accessories**  
 3001451 Automix Dual-barrel Syringe Tips (25)

**premier®** Inspired solutions for daily dentistry®



Premier® Dental Products Company • 1710 Romano Drive, Plymouth Meeting, PA 19462, USA  
 Phone 888-670-6100 • 610-239-6000 • Fax 610-239-6171 • premierdentalco.com • E-mail: info@premierdentalco.com

Lit# 8621310 1124026 Rev0 NKP5M